


HADEX: AN R PACKAGE AND WEB-SERVER FOR ANALYSIS OF DATA FROM HDX-MS EXPERIMENTS


DOMINIK CYSEWSKI

INSTITUTE OF BIOCHEMISTRY AND BIOPHYSICS
POLISH ACADEMY OF SCIENCES




The longer incubation time, the more protected hydrogens are being replaced by deuters.

APPLICATIONS



G. R. Masson et. al. (2019). Recommendations for performing, interpreting and reporting hydrogen deuterium exchange mass spectrometry (HDX-MS) experiments. *Nature Methods*, 16(7):595–602..


MULTI-STATE ANALYSIS



Peptides may come from proteins in different states, i.e. bounded by different cofactors.

During the analysis of the HDX-MS data we have to consider both time of the incubation and position of a peptide in a sequence.

RELATIVE DEUTERATION



DEUTERATION KINETICS

Start Woods plot Coverage Sequence data **Kinetics** Summary Report



Select parameters for the plot.

☐ Theoretical calculations

Choose values type:

- ☒ relative
☐ absolute

Choose time parameters:


IN OUT
0.001 1440

Choose peptide:

Sequence	State	Start	End
INITSSASQECTRLN	CD160	1	15
INITSSASQECTRLN	CD160_HVEM	1	15
LICTVW	CD160	16	21
LICTVW	CD160_HVEM	16	21
LICTVWHKKEEAEG	CD160	16	29
LICTVWHKKEEAEG	CD160_HVEM	16	29
ICTVWHKKEEAEG	CD160	17	29
ICTVWHKKEEAEG	CD160_HVEM	17	29
FVVFLL	CD160	30	34
FVVFLL	CD160_HVEM	30	34

Kinetic plot

Data



Save chart (.svg)

REPORTING

Start

Woods plot


Coverage

Sequence data

Kinetics

Summary

Report



Name	Value
HDX time course	8
Number of peptides	41
Sequence coverage	0.8636
Average peptide length	14.6829
Redundancy	4.5606
Replicates	4
Significant differences in HDX	0.98 0.99

All

All

Excel

PDF

Showing 1 to 7 of 7 entries.

This table was adapted from the HDX Data Summary (Mason, G.R., Burke, J.E., Afari, N.G., Anand, C.S., Borchers, C., Drier, S., Dou-Assaf, G.M., Engen, J.R., Engländer, S.W., Folmer, J., et al. (2019). Recommendations for performing, interpreting and reporting hydrogen deuterium exchange mass spectrometry (HDX-MS) experiments. Nat Methods 16, 595-602).

Start

Woods plot


Coverage

Sequence data

Kinetics

Summary

Report



Choose items for report:

☒ Position Frequency

☒ Peptide Coverage

☒ Comparison Plot

☒ Theoretical Comparison Plot

☒ Woods Plot

☒ Theoretical Woods Plot

☒ Kinetic Plot

☐ Theoretical Kinetic Plot

☐ Position Frequency Data

☐ Peptide Coverage Data

☐ Comparison Plot Data

☐ Theoretical Comparison Plot Data

☐ Woods Plot Data

☐ Theoretical Woods Plot Data

☐ Kinetic Plot Data

☐ Theoretical Kinetic Plot Data

Elements chosen for report have the same parameters as chosen in panel e.g. axis range and title. Adjust parameters for plots as needed in the report.

Create report!

SUMMARY AND AVAILABILITY

Summary:

- rapidly developing technologies require a flexible framework,
- the methods of data analysis should follow the development of both technology and expectations of its users.

Availability:

- a web-server (<http://mslab-ibb.pl/shiny/HaDeX/>),
- the R package (<https://CRAN.R-project.org/package=HaDeX>),
- a standalone software (<https://sourceforge.net/projects/HaDeX/>).

ACKNOWLEDGEMENTS

- Mass Spectrometry Lab, Institute of Biochemistry and Biophysics, PAS.
- MI² Data Lab, Faculty of Mathematics and Information Science, Warsaw University of Technology.



Funding: Foundation of Polish Science TEAM TECH CORE FACILITY/2016-2/2 Mass Spectrometry of Biopharmaceuticals - improved methodologies for qualitative, quantitative and structural characterization of drugs, proteinaceous drug targets and diagnostic molecules.

ACKNOWLEDGEMENTS

HaDeX developers:

- Weronika Puchala
(main developer),
- Dominik Rafacz
(frontend developer),
- Michał Burdukiewicz.

